



MODU-MAX CTS

sterile instrument storage

PRODUCT SPOTLIGHT

FIRST-IN FIRST-OUT FIFO Label

✓ Highly Visible

✓ Set In Frame

✓ Easy Rotation

How FIFO Labels Work

The **FIRST-IN FIRST-OUT (FIFO)** label mounts directly to MODU-MAX CTS frames. MODU-MAX CTS with FIFO labels keeps your SPD compliant, organized, and safe.

1 Green Side Signals 'Use First'

Flip the bi-color disc to signal the **oldest set should be received first**.


When **Green** side faces outward, staff knows this shelf holds the oldest sterilized equipment-**Retrieve these first**.



2 Retrieve & Rotate to Red

After retrieving the blue wrap set, rotate disc back to **Red**.

Red is a clear visual cue that eliminates guesswork and upholds sterility standards



3 AAMI ST70 Compliant

Proper rotation **reduces risk of expired blue wrap sets** reaching the OR.



STOP STACKING



START COMPLYING

AAMI ST79 STORAGE GUIDELINES

Why Storage with FIFO Labels Matter

11.1.1

NO STACKING

Stacking can result in damage to the wrap caused by undue pressure from the weight. 11.1.1

11.1.1

NO COMPRESSION

Compression of packages can force air and microorganisms into the package...which lead to contamination. 11.1.1

11.2.1

HANDLE WITH CARE

Sterile items should be handled to avoid: **sliding, dragging, crushing, bending, compressing, puncturing...** compromising the sterility of the contents. 11.2.1

ANSI/AAMI ST79:2017 American National Standard: Comprehensive guide to steam sterilization assurance in health care facilities Published by AAMI 4301 N. Fairfax Dr., Suite 301 Arlington, VA 22203-1633 www.aami.org © 2017 by the Association for the Advancement of Medical Instrumentation